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Available online <http://icr.toray.com/content/11/2/R47>

Research

Effectiveness of polymyxin B-immobilized fiber column in sepsis: a systematic review

Dima N Chai¹, Mark A Penzella², Rinaldo Bultroni³, Massimo de Cal⁴, Natalia Polanco¹, Valentina Corradi¹, Paolo Lentini¹, Federico Nalesso¹, Takuya Ueno⁵, V Marco Ranieri⁶ and Claudio Ronco¹¹Department of Nephrology, Ospedale San Bortolo, Viale Rodari 37, 36100 Verona, Italy; ²Section of Nephrology, Department of Medicine, Yale University School of Medicine, 333 Cedar Street PMP 137, New Haven, CT 06520, USA; ³Department of Internal Care and Department of Medicine, Alfred & Redfern Medical Centre, Studley Road, Mordialloc, Victoria 3008, Australia; ⁴Transplantation Unit, Surgical Services, Massachusetts General Hospital, 55 Fruit Street White St, Boston, MA 02114, USA; ⁵Department of Anesthesia and Intensive Care, Ospedale San Giovanni-Battista, Corso Bramante 18, 10128 Torino, ItalyCorresponding author: Claudio Ronco, ronco@ospedale

Received: 4 Jan 2007; Revisions requested: 9 Feb 2007; Revisions received: 1 Mar 2007; Accepted: 20 Apr 2007; Published: 29 Apr 2007

Critical Care 2007, 11:R47 (doi:10.1186/cc5780)

This article is online at <http://icr.toray.com/content/11/2/R47>

© 2007 Chai et al; licensee BioMed Central Ltd.

This is an open access article distributed under the terms of the Creative Commons Attribution License (<http://creativecommons.org/licenses/by/2.0>), which permits unrestricted use, distribution, and reproduction in any medium, provided the original work is properly cited.

Abstract

Introduction Severe sepsis and septic shock are common problems in the intensive care unit and carry a high mortality. Endotoxin, one of the principal components on the outer membrane of gram-negative bacteria, is considered important to their pathogenesis. Polymyxin B bound and immobilized to polystyrene fibers (PMX-F) is a medical device that aims to remove circulating endotoxin by adsorption, theoretically preventing the progression of the biological cascade of sepsis. We performed a systematic review to describe the effect in septic patients of direct hemoperfusion with PMX-F on outcomes of blood pressure, use of vasopressor drugs, oxygenation, and mortality reported in published studies.

Methods We searched PubMed, the Cochrane Collaboration Database, and bibliographies of retrieved articles and consulted with experts to identify relevant studies. Prospective and retrospective observational studies, pre- and post-intervention design, and randomized controlled trials were included. Three authors reviewed all citations. We identified a total of 28 publications - 8 randomized controlled trials, 7 non-randomized parallel studies, and 12 pre-post design studies - that reported at least one of the specified outcome measures: pooled sample

size, 1,425 patients; 978 PMX-F and 447 conventional medical therapy).

Results Overall, mean arterial pressure (MAP) increased by 19 mm Hg (95% confidence interval [CI], 15 to 22 mm Hg, $p < 0.001$), representing a 26% mean increase in MAP (range, 14% to 42%), whereas dopamine/dobutamine dose decreased by 1.8 µg/kg per minute (95% CI, 0.4 to 3.3 µg/kg per minute, $p = 0.01$) after PMX-F. There was significant intertrial heterogeneity for these outcomes ($p < 0.001$), which became non-significant when analysis was stratified for baseline MAP. The mean arterial partial pressure of oxygenation of inspired oxygen (P_{IO}/F_O₂) ratio increased by 32 units (95% CI, 23 to 41 units, $p < 0.001$). PMX-F therapy was associated with significantly lower mortality risk (risk ratio, 0.53; 95% CI, 0.43 to 0.65). The trials assessed had suboptimal method quality.

Conclusion Based on this critical review of the published literature, direct hemoperfusion with PMX-F appears to have favorable effects on MAP, dopamine use, P_{IO}/F_O₂ ratio, and mortality. However, publication bias and lack of blinding need to be considered. These findings support the need for further rigorous study of this therapy.

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
Postoperative Polymyxin B Hemoperfusion and Mortality in Patients With Abdominal Septic Shock: A Propensity-Matched Analysis*

Masao Iwagami, MD, MPH¹; Hideo Yasunaga, MD, PhD²; Kent Doi, MD, PhD³; Hiromasa Horiguchi, PhD⁴; Kiyohide Fushimi, MD, PhD⁵; Takehiro Matsubara, MD, PhD⁶; Naoki Yahagi, MD, PhD⁷; Eisei Noiri, MD, PhD¹

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Case Report

Hemoperfusion Treatment in a Septic Shock Patient with Autosomal Dominant Polycystic Kidney Disease and Increased HMGB1 Protein Levels

Takana Nakamura¹, Naohito Kawaguchi², Norihiko Ueda³, Shingo Yamada⁴, Hiroyuki Kikuchi⁵

¹Department of Medicine, Shimizu City General Hospital, Shimizu, Japan; ²Department of Pathology, Shimizu City General Hospital, Shimizu, Japan; ³Department of Radiology, Shimizu City General Hospital, Shimizu, Japan; ⁴Department of Anesthesiology, Shimizu City General Hospital, Shimizu, Japan; ⁵Department of Laboratory Medicine, Shimizu City General Hospital, Shimizu, Japan

Key Words: Hemoperfusion, septic shock, autosomal dominant polycystic kidney disease, HMGB1, protein levels

Abstract: This case report describes a patient with septic shock caused by peritonitis. The patient was treated with hemoperfusion (HP) and received a large volume of plasma infusion. The patient was discharged on day 23 after admission.

Introduction: Septic shock and severe organ dysfunction are common problems in the intensive care unit, and the associated mortality rate is high. (1) In 1991, a patient was administered HP (PMX-1) column was supported by the Japanese Society for Intensive Care Medicine (JICM). (2) In a recent study, HP (PMX-1) treatment was shown to be effective in septic shock patients. (3) However, the effect of PMX-1 treatment in septic shock patients with a blood evidence level of more than 1,000 pg/mL is unclear. High-mobility group box 1 (HMGB1) is a secreted DNA-binding protein that acts as a proinflammatory cytokine and inhibits organ function.

Case Report: A 68-year-old female patient with a long history of autosomal dominant polycystic kidney disease (ADPKD) was admitted to the intensive care unit (ICU) with septic shock caused by peritonitis. The patient had a blood evidence level of HMGB1 of 2,200 pg/mL. She was treated with HP (PMX-1) and received a large volume of plasma infusion. The patient was discharged on day 23 after admission.

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