

表 397 ETRF装着率(%) ET測定値(EU/mL)別 (ベッドサイドコンソール台数≥1)

ET測定値	0(装着せず)	10未満	10～	20～	30～	40～	50～	60～	70～	80～	90～	100(全台装着)	合計	記載なし	総計	平均	標準偏差
0.001EU/mL未満	42	37	34	35	39	22	35	39	45	72	142	2,623	3,165	2	3,167	92.75	21.39
(%)	(1.3)	(1.2)	(1.1)	(1.1)	(1.2)	(0.7)	(1.1)	(1.2)	(1.4)	(2.3)	(4.5)	(82.9)	(100.0)				
0.001～0.01EU/mL未満	30	20	19	16	15	7	7	9	8	9	25	381	546	2	548	81.19	34.41
(%)	(5.5)	(3.7)	(3.5)	(2.9)	(2.7)	(1.3)	(1.3)	(1.6)	(1.5)	(1.6)	(4.6)	(69.8)	(100.0)				
0.01～0.05EU/mL未満	18	7	9	6	5	3	5	3	6	7	11	131	211		211	77.17	36.71
(%)	(8.5)	(3.3)	(4.3)	(2.8)	(2.4)	(1.4)	(2.4)	(1.4)	(2.8)	(3.3)	(5.2)	(62.1)	(100.0)				
0.05～0.1EU/mL未満	8	4	1	1		2	2	1	3	3	2	36	63		63	73.25	39.36
(%)	(12.7)	(6.3)	(1.6)	(1.6)		(3.2)	(3.2)	(1.6)	(4.8)	(4.8)	(3.2)	(57.1)	(100.0)				
0.1～0.25EU/mL未満	4	2	2	2	1	1	1			1		34	48		48	77.00	38.34
(%)	(8.3)	(4.2)	(4.2)	(4.2)	(2.1)	(2.1)	(2.1)			(2.1)		(70.8)	(100.0)				
0.25～0.5EU/mL未満	1			3			2	1	1	1		11	20		20	75.20	33.39
(%)	(5.0)			(15.0)			(10.0)	(5.0)	(5.0)	(5.0)		(55.0)	(100.0)				
0.5EU/mL以上			1			2	1	1		2	2	13	22		22	85.29	25.17
(%)			(4.5)			(9.1)	(4.5)	(4.5)		(9.1)	(9.1)	(59.1)	(100.0)				
合計	103	70	66	63	60	37	53	54	63	95	182	3,229	4,075	4	4,079	89.78	25.77
(%)	(2.5)	(1.7)	(1.6)	(1.5)	(1.5)	(0.9)	(1.3)	(1.3)	(1.5)	(2.3)	(4.5)	(79.2)	(100.0)				
不明	40	3	5	6	1		9	6	5	2	5	81	163	1	164	62.73	43.65
(%)	(24.5)	(1.8)	(3.1)	(3.7)	(0.6)		(5.5)	(3.7)	(3.1)	(1.2)	(3.1)	(49.7)	(100.0)				
記載なし	16		1	2	1		1	2	1		1	29	54	7	61	62.01	45.63
(%)	(29.6)		(1.9)	(3.7)	(1.9)		(1.9)	(3.7)	(1.9)		(1.9)	(53.7)	(100.0)				
総計	159	73	72	71	62	37	63	62	69	97	188	3,339	4,292	12	4,304	88.41	27.64
(%)	(3.7)	(1.7)	(1.7)	(1.7)	(1.4)	(0.9)	(1.5)	(1.4)	(1.6)	(2.3)	(4.4)	(77.8)	(100.0)				

数値下のかっこ内は行方向の合計に対する%です。

(C)Japanese Society for Dialysis Therapy

2014年調査