

表 481 ET測定値(EU/mL) 細菌数測定の高度別 (ベッドサイドコンソール台数 1)

細菌数測定の高度	0.001 EU/mL未満	0.001～0.01 EU/mL未満	0.01～0.05 EU/mL未満	0.05～0.1 EU/mL未満	0.1～0.25 EU/mL未満	0.25～0.5 EU/mL未満	0.5 EU/mL以上	合計	不明	記載なし	総計
なし	55	17	10		5		1	88	53	43	184
(%)	(62.5)	(19.3)	(11.4)		(5.7)		(1.1)	(100.0)			
毎日	11	2						13			13
(%)	(84.6)	(15.4)						(100.0)			
毎週	120	17	4	4	1	2		148			148
(%)	(81.1)	(11.5)	(2.7)	(2.7)	(0.7)	(1.4)		(100.0)			
隔週	197	15	10	2	1	2	2	229			229
(%)	(86.0)	(6.6)	(4.4)	(0.9)	(0.4)	(0.9)	(0.9)	(100.0)			
毎月	2,399	326	101	24	17	8	8	2,883	3		2,886
(%)	(83.2)	(11.3)	(3.5)	(0.8)	(0.6)	(0.3)	(0.3)	(100.0)			
年数回	263	53	15	7	4	5	1	348	6		354
(%)	(75.6)	(15.2)	(4.3)	(2.0)	(1.1)	(1.4)	(0.3)	(100.0)			
年1回	355	26	17	8	8	3	4	421	4		425
(%)	(84.3)	(6.2)	(4.0)	(1.9)	(1.9)	(0.7)	(1.0)	(100.0)			
合計	3,400	456	157	45	36	20	16	4,130	66	43	4,239
(%)	(82.3)	(11.0)	(3.8)	(1.1)	(0.9)	(0.5)	(0.4)	(100.0)			
不明	6	3	3	1	1			14	62	2	78
(%)	(42.9)	(21.4)	(21.4)	(7.1)	(7.1)			(100.0)			
記載なし										1	1
(%)											
総計	3,406	459	160	46	37	20	16	4,144	128	46	4,318
(%)	(82.2)	(11.1)	(3.9)	(1.1)	(0.9)	(0.5)	(0.4)	(100.0)			

数値下のかっこ内は行方向の合計に対する%です。

##C0055##

(C)Japanese Society for Dialysis Therapy

2016年調査